

**Table 1****Sensitivity of AC2 compared with routine culture in identifying *Neisseria gonorrhoeae***

|                               | Sensitivity |        | Culture (%) | 95% CI |
|-------------------------------|-------------|--------|-------------|--------|
|                               | AC2 (%)     | 95% CI |             |        |
| Men                           |             |        |             |        |
| With urethral symptoms        | 99          | 95-100 | 79          | 71-85  |
| Without urethral symptoms     | 94          | 69-100 | 29          | 11-56  |
| With rectal symptoms          | 100         | 66-100 | 60          | 27-86  |
| Without rectal symptoms       | 92          | 62-100 | 54          | 26-80  |
| Women                         |             |        |             |        |
| With vulvovaginal symptoms    | 100         | 95-100 | 53          | 38-68  |
| Without vulvovaginal symptoms | 100         | 87-100 | 47          | 30-65  |